



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Simulator, childbirth, complete, dark

**Product Code :**  
LBNY-0015-220015



**Description :**

Complete childbirth simulator, dark complexion, used for demonstration and practice of normal childbirth of single foetus and early postpartum management.

Simulation will cover the following obstetric procedures:

- Palpation of foetal fontanelles;
- Normal vaginal delivery of single foetus;
- Cutting and clamping umbilical cord;
- Examination of placenta, umbilical cord and membranes during normal delivery
- Early postpartum management;

**Technical Specification :**

Environmental conditions:

Storage conditions:

Storage humidity max 95% RH

Storage temperatures -10 to 50 °C (14 - 122 °F)

Operating conditions:

Operating humidity 5 - 100% RH

Operating temperatures 15 to 40 °C (59 - 104 °F)

---

# Laboratory Instrument India

