



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Simulator, childbirth, complete, light

**Product Code :**  
LBNY-0015-2200016



**Description :**

Complete childbirth simulator, light complexion, used for demonstration and practice of normal childbirth of single foetus and early postpartum management.

Normal vaginal delivery of single foetus;

Simulation will cover the following obstetric procedures:

Examination of placenta, umbilical cord and membranes during normal delivery

Cutting and clamping umbilical cord;

Early postpartum management;

Palpation of foetal fontanel;

**Technical Specification :**

Environmental considerations:

Storage conditions:

Storage humidity max 95% RH

Storage temperatures -10 to 50 °C (14 - 122 °F)

Operating conditions:

Operating humidity 5 - 100% RH

Operating temperatures 15 to 40 °C (59 - 104 °F)

Supplied with:

Removable palpable abdominal cover and moveable parts for positioning foetal baby/babies;

1 umbilical cord, that can be cut and attached several times for cordcutting and clamping demonstration

2 umbilical ties

---

One foetus as specified in general design requirements;  
1 placenta with membranes;  
Instruction for use in English, Spanish and French, including pictographs;  
A soft/hard carrying case with durability guaranteed under normal use.

Weight and volume:

Volume per unit: 0.058 cbm

Unit qty in one carton: 1 piece

Gross weight per carton: 5.4 kg

Volume per carton: 0.058 cbm

Carton size: 0.63 x 0.33 x 0.28 m

Gross weight per unit: 5.4 kg.



**Laboratory Instrument India**