



Email : sales@laboratoryinstrumentindia.com

Product Name :
Simulator, newborn complete care, light

Product Code :
LBNY-0015-2200018



Description :

Newborn simulator to facilitate learning complete newborn care, including cord cutting and ligation, as well as newborn resuscitation, light complexion
Simulator consists of a light skin manikin, that has an appearance and size of full term infant with umbilical cord.
Material of the simulator: Thermoplastic polyurethane - TPU (primary material for body skin), Polyvinylchloride - PVC (primary material in chest and face skin) silicone rubber, Polyoxymethylene - POM, Polypropylene - PP, and iron/steel.
Texture of the simulator: soft and smooth, friction free to demonstrate the desired procedures. Realistically sculpted to resemble and simulate a baby, anatomically accurate and feeling smooth / resilient / bony as relevant and suitable for simulation.

Technical Specification :

Supplied with:
1 umbilical cord, that can be cut and attached several times for demonstration purposes
2 sheets for drying and wrapping a newborn.
1 head cap
2 umbilical ties
Suction device, portable, hand-held, manual to clear newborn's airway.
Shoulder bag for transport/storage
Translucent, made of silicone or other similar material, which is easy to clean and can withstand sterilization
Self-inflating neonatal resuscitation bag with 2 translucent masks: size 0 and size 1.

Environmental conditions:

Storage conditions with water: +20°C / 95% RH for 3 years. (Short term for 50°C - 2 weeks)

Operating conditions:

Operating humidity: 5 to 95% RH

Operating temperature: +15°C to +40°C

Weight and volume:

Volume per unit: 0.006 cbm

Gross weight per unit: 1.2 kg.



Laboratory Instrument India