



Email : sales@laboratoryinstrumentindia.com

Product Name :
Trainer, breastfeeding, dark

Product Code :
LBNY-0015-2200021



Description :

Breastfeeding trainer is used to support demonstration of proper positioning and attachment during breastfeeding and learning breast milk expression.

Trainer is attachable to manikin or individual demonstrating the procedures.

Materials of the simulator: Neoprene, poly vinyl chloride (PVC),

Polypropylene - PP, Silicone, Chlorinated polyvinyl chloride- -CPVC, Polyester, Nylon, Lycra,

Wearable task trainer to simulate adult female breasts, has a realistic appearance and size of nursing mother's breasts.

Texture of the simulator: soft and smooth, friction free to demonstrate the desired procedures. Realistically sculpted and anatomically accurate, feeling smooth / resilient, as relevant and malleable parts suitable for simulation.

Technical Specification :

Supplied with:

Attachments to ensure the model can be secured on a mannequin or demonstrator;

A soft/hard carrying case with durability guaranteed under normal use.

Premature infant feeding cup (paladai).

Instruction for use in English, Spanish and French with pictographs

Appropriate refillable tanks.

Environmental conditions:

Storage conditions:

Storage environmental conditions: +20°C for 3 years.

(Short term for 50°C for up to 2 weeks)

Operating humidity: 5 to 95% RH

Operating conditions: Operating temperature: +15°C to +40°C

Dimensions:

Volume per unit(M3): 0.004

Gross weight per unit: 0.548 kg.



Laboratory Instrument India