




Email : sales@laboratoryinstrumentindia.com

Product Name :
Receptacle,waste,pedal action,20L

Product Code :
LBNY-0016-220001



Description :

Round waste receptacle with foot pedal opens self-closing lid, stainless steel.
With removable inner heavy gauge plastic bin with handle.
Receptacle, waste, pedal action, 20 litres
Capacity: 19-20 litres.
Height: 38-45 cm.

Technical Specification :

Weight/Volume/Dimensions:
Estimated weight: 6.8 kg per unit
Estimated volume: 0.192 m3 per carton



Laboratory Instrument India