



Email : sales@laboratoryinstrumentindia.com

Product Name :
Table Vice 1000

Product Code :
LBNY-0004-820018



Description :

Table Vice 1000

Technical Specification :

Table Vice 1000

Made from grey iron castings - a very strong material with excellent shock absorbing capacity and immense strength under compression, these are machined to close tolerances, are designed for smooth operation and reliability.

The front jaw beam fully encloses the vice spindle protecting the threads from contaminants and damage.

Spindle is forged with square threads for better clamping strength.

The floating; unbreakable s.g iron nut securely locked in the body ensures thread surface and constant alignment between body and the sliding jaw.

Vice jaws are hardened for increased resistance to wear; are knurled for better grip.

Vices have a powder coated finish with main screw and handle plated.

4-1/2-Inch swivel base Vice

4.5" Jaw Width

3" Jaw Opening

Holds 0.6" D to 1.85" D pipes and tubes

240-Degree Swivel Base.

Laboratory Instrument India

