



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Basin,kidney,polypropylene,475ml

**Product Code :**  
LBNY-0016-220002



**Description :**

Material: Polypropylene, autoclavable.  
Capacity: approx. 475ml  
Length: approx. 205mm  
Width: approx. 95mm  
Height: approx. 45mm.

**Technical Specification :**

Weight/Volume/Dimensions:  
Estimated volume: 0.264cdm  
Estimated weight: 0.054kg.



**Laboratory Instrument India**