



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Basin,kidney,stainless steel,825ml

**Product Code :**

LBNY-0016-220003



**Description :**

Basic container for soiled dressings during nursing and surgical care.

Smooth surface.

Material: Austenitic stainless steel.

Austenitic stainless steel composition: 18 - 20% chromium, 8 - 10 % nickel.

Reusable.

Container, kidney shaped.

Capacity: approx. 800 to 850ml.

Thickness: approx. 0.75-0.85mm.

Length: approx. 240-260 mm.

Width: approx. 130-150mm.

Height: approx. 30-50mm.

**Technical Specification :**

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.

---

# Laboratory Instrument India

