




Email : sales@laboratoryinstrumentindia.com

Product Name :

Bowl,round,polypropylene,6 - 8L

Product Code :

LBNY-0016-220004



Description :

Round bowl.

Material: Polypropylene, autoclavable.

Capacity: approximately 6 - 8 litres.

Diameter: approximately 346mm.

Height: approximately 120mm.

Technical Specification :

Weight/Volume/Dimensions:

Estimated volume: 2.266cdm

Estimated weight: 0.300kg.



Laboratory Instrument India