




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Bowl,round,stainless steel,4L

**Product Code :**  
LBNY-0016-220005



**Description :**

Container, round.  
Smooth surface.  
Material: Austenitic stainless steel.  
Austenitic stainless steel composition: 18 - 20% chromium, 8 - 10% nickel.  
Reusable.  
Capacity: approx. 4 litres.  
Internal diameter: approx. 300 mm.  
Height: approx. 100 mm.  
Thickness: approx. 0.8 mm.

**Technical Specification :**

Environmental conditions:  
Humidity: 30% RH to 75% RH.  
Storage and transport conditions: 1°C to 40°C.



**Laboratory Instrument India**

