



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bowl, stainless steel, 600ml

Product Code :
LBNY-0016-220007



Description :

Container, round.
Smooth surface.
Material: Austenitic stainless steel.
Austenitic stainless steel composition: 18-20% chromium, 8-10% nickel.
Reusable.
Capacity: approx. 550 - 650 ml.
Inner diameter: approx. 115-135 mm.
Height: approx. 50-70 mm.
Thickness: approx. 0.75-0.85 mm.

Technical Specification :

Environmental conditions:
Humidity: 30% RH to 70% RH.
Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

