



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Bowl, stainless steel, 600ml

**Product Code :**  
LBNY-0016-220007



**Description :**

Container, round.  
Smooth surface.  
Material: Austenitic stainless steel.  
Austenitic stainless steel composition: 18-20% chromium, 8-10% nickel.  
Reusable.  
Capacity: approx. 550 - 650 ml.  
Inner diameter: approx. 115-135 mm.  
Height: approx. 50-70 mm.  
Thickness: approx. 0.75-0.85 mm.

**Technical Specification :**

Environmental conditions:  
Humidity: 30% RH to 70% RH.  
Storage and transport conditions: 1°C to 30°C.



**Laboratory Instrument India**

