



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bowl, stainless steel, 180ml

Product Code :
LBNY-0016-220008



www.laboratoryinstrumentindia.com

Description :

Basic container for nursing and surgical care.
Material: Austenitic stainless steel.
Austenitic stainless steel composition: 18-20% chromium, 8-10% nickel.
Reusable.
Container, round.
Smooth surface.
Capacity: approx. 170 to 190ml.
Inner diameter: approx. 75-85 mm.
Height: approx. 36-48 mm.
Thickness: approx. 0.75-0.85mm.

Technical Specification :

Environmental conditions:
Humidity: 30% RH to 75% RH.
Storage and transport conditions: 1°C to 40°C.



Laboratory Instrument India

