



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Tray,instr,ss,225x125x50mm,w/cover

**Product Code :**

LBNY-0016-220009



**Description :**

Tray, instruments.

Material: Austenitic stainless steel.

Seamless tray with cover.

Rectangular with rounded corners, smooth surface.

Austenitic stainless steel composition: 18-20% chromium, 8-10% nickel.

Reusable.

Length: approx. 215-235 mm.

Width: approx. 120-130 mm.

Height: approx. 45-50 mm.

Thickness: approx. 0.75-0.85 mm.

**Technical Specification :**

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.



**Laboratory Instrument India**

