



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tray,instr,ss,310x195x60mm,w/cover

Product Code :

LBNY-0016-2200010



Description :

Tray, instruments.

Material: Austenitic stainless steel.

Austenitic stainless steel composition: 18-20% chromium, 8-10% nickel.

Reusable.

Seamless tray with cover.

Rectangular with rounded corners, smooth surface.

Length: approx. 300-320.mm.

Width: approx. 185-205.mm.

Height: approx. 55-65.mm.

Thickness: approx. 0.75-0.85.mm.

Technical Specification :

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

