



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Tray,dressing,ss,480x330x20mm

**Product Code :**

LBNY-0016-2200011



**Description :**

Tray, dressing.

Material: Austenitic stainless steel.

Seamless tray, rectangular with rounded corners, smooth surface.

Reusable.

Length: approx. 460-500 mm.

Width: approx. 315-345 mm.

Height: approx. 19-21 mm.

Thickness: approx. 0.75-0.85 mm.

Austenitic stainless steel composition: 18-20% chromium, 8-10% nickel.

**Technical Specification :**

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.



**Laboratory Instrument India**

