



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tray,dressing,ss,480x330x20mm

Product Code :

LBNY-0016-2200011



Description :

Tray, dressing.

Material: Austenitic stainless steel.

Seamless tray, rectangular with rounded corners, smooth surface.

Reusable.

Length: approx. 460-500 mm.

Width: approx. 315-345 mm.

Height: approx. 19-21 mm.

Thickness: approx. 0.75-0.85 mm.

Austenitic stainless steel composition: 18-20% chromium, 8-10% nickel.

Technical Specification :

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

