



Email : sales@laboratoryinstrumentindia.com

Product Name :
Scissors Mechanical 100

Product Code :
LBNY-0004-820019



Description :

Scissors Mechanical 100

Technical Specification :

Scissors Mechanical 100

This scissors head suitable for cut off metal sheet, wire mesh, aluminium plate, stainless steel plate, plastic, cloth, leather etc.

Straight cut airline scissors

Material; chrome vanadium forged steel

Cuttingrange: 18GA-22GA

Cutable thickness: Cold rolled steel: 18GA /

Stainless steel: 22GA

Blade Length; 42mm



Laboratory Instrument India