



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Tray,dressing,ss,300x200x30mm

**Product Code :**

LBNY-0016-2200012



**Description :**

Tray, dressing, stainless steel.

Seamless tray, rectangular with rounded corners, smooth surface.

Material: Austenitic stainless steel.

Length: approx. 285-315 mm.

Width: approx. 190-220 mm.

Height: approx. 28-32 mm.

Thickness: approx. 0.75-0.85 mm.

Austenitic stainless steel composition: 18-20% chromium, 8-10% nickel.

Reusable.

**Technical Specification :**

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.



**Laboratory Instrument India**

