



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tray, dressing, ss, 300x200x30mm

Product Code :

LBNY-0016-2200012



Description :

Tray, dressing, stainless steel.
Seamless tray, rectangular with rounded corners, smooth surface.
Material: Austenitic stainless steel.
Length: approx. 285-315 mm.
Width: approx. 190-220 mm.
Height: approx. 28-32 mm.
Thickness: approx. 0.75-0.85 mm.
Austenitic stainless steel composition: 18-20% chromium, 8-10% nickel.
Reusable.

Technical Specification :

Environmental conditions:
Humidity: 30% RH to 70% RH.
Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

