



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Jar,thermometer,pp,11cm,w/o cover

**Product Code :**

LBNY-0016-2200014



**Description :**

Material: Polypropylene, autoclavable.  
Overall height: approximately 120mm.  
Depth: approximately 110mm.  
Internal diameter: approximately 30mm.

**Technical Specification :**

Weight/Volume/Dimensions:  
Estimated volume: 1.573cdm  
Estimated weight: 0.475kg



**Laboratory Instrument India**