



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Bottle,plastic,250ml,wash bottle

**Product Code :**

LBNY-0016-2200018



**Description :**

Container: Wash bottle, round, mouth with smaller diameter than body, with screw closure and riser tube. Bottle assembled.

Bottle, plastic, 250ml, wash bottle suitable that will not degrade with long term storage of chlorine and iodine solutions in hospital environment.

Material: Plastic.

Capacity: 250 ml.

**Technical Specification :**

Weight and volume:

Estimated volume: 1.075cdm

Estimated weight: 0.094kg



**Laboratory Instrument India**