



Email : sales@laboratoryinstrumentindia.com

Product Name :

Prongs,nasal,Oxygen,adult,s.u.

Product Code :

LBNY-0017-2200013



www.laboratoryinstrumentindia.com

Description :

Nasal cannula (nasal prongs), device designed for easy administration of oxygen into the patient nose through two small prongs placed in the nostrils, providing comfort for the patient.

The device consists of:

Star lumen main tube to avoid accidental blockage.

Over-the-Ear tubing adjustable

Soft funnel shaped connector to facilitate easy connection to oxygen source.

Adjustable, smoothly finished nasal tips for maximum patient comfort.

Soft twin prongs nasal tips to ensure equal oxygen flow to both.

Length of the oxygen tube: +/-2 m.

Material: Preferably soft and kink resistant polyvinyl chloride (PVC).

Single use.

Non-sterile.

Size selected: Adult

Technical Specification :

WEIGHT/VOLUME (packaged):

Estimated volume: 0.53cdm

Estimated weight: 0.091kg



Laboratory Instrument India