



Email : sales@laboratoryinstrumentindia.com

Product Name :

Prongs,nasal,Oxygen,child,s.u.

Product Code :

LBNY-0017-2200014



www.laboratoryinstrumentindia.com

Description :

Nasal cannula (nasal prongs), device designed for easy administration of oxygen into the patient nose through two small prongs placed in the nostrils, providing comfort for the patient.

The device consists of:

Material: Preferably soft and kink resistant polyvinyl chloride (PVC).

Star lumen main tube to avoid accidental blockage.

Over-the-Ear tubing adjustable

Soft funnel shaped connector to facilitate easy connection to oxygen source.

Soft twin prongs nasal tips to ensure equal oxygen flow to both.

Adjustable, smoothly finished nasal tips for maximum patient comfort.

Length of the oxygen tube: +/- 2 m.

Size selected: Child

Single use.

Non-sterile.

Technical Specification :

WEIGHT/VOLUME (packaged):

Estimated volume: 0.53cdm

Estimated weight: 0.091kg.

Laboratory Instrument India

