



Email : sales@laboratoryinstrumentindia.com

Product Name :

Prongs,nasal,Oxygen,neonate,s.u.

Product Code :

LBNY-0017-220015



www.laboratoryinstrumentindia.com

Description :

Nasal cannula (nasal prongs), device designed for easy administration of oxygen into the patient nose through two small prongs placed in the nostrils, providing comfort for the patient.

The device consists of:

Soft twin prongs nasal tips to ensure equal oxygen flow to both.

Soft funnel shaped connector to facilitate easy connection to oxygen source.

Length of the oxygen tube: +/- 2 m.

Material: Preferably soft and kink resistant polyvinyl chloride (PVC).

Adjustable, smoothly finished nasal tips for maximum patient comfort.

Star lumen main tube to avoid accidental blockage.

Over-the-Ear tubing adjustable

Size selected: Neonate

Non-sterile.

Single use.

Technical Specification :

WEIGHT/VOLUME (packaged):

Estimated volume: 0.53 cdm

Estimated weight: 0.091 kg

Laboratory Instrument India

