



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Pump,suction,foot-operated

**Product Code :**  
LBNY-0017-2200016



www.laboratoryinstrumentindia.com

#### **Description :**

Basic hospital equipment for health structures and emergency situations in wards, emergency room, operating theatre, delivery room, intensive care unit, etc.

The suction pump should be operated only by a person who has received adequate training in pharyngeal and tracheal suction techniques.

High performance suction pump, foot-operated for pharyngeal and tracheal suction and during minor surgical interventions. Supplied with instruction manual and diagrams covering the function of the pump, how to use it, dismantle and assemble it, to clean, disinfect and sterilize it, its maintenance and spare-parts.

#### **Technical Specification :**

The unit provides suction a foot operated pedal.

The unit is provided with a suction tube of at least 1.3 meters.

The suction tube is made of medical grade silicone which is autoclavable at 121°C.

The user adjustable suction ranges between 0 and at least 525 mmHg (70kPa).

The flow rate at maximum vacuum is 40 L/min (2.4 m<sup>3</sup>/h).

The unit is mounted on a stable platform and easy to carry.

The base has an anti-slip provision, preventing the pump from sliding away when in use.

The unit has maintenance free operating requiring little or no oil.

The pump is equipped with a provision preventing fluid flowing to the pump if the collection jar is full.

The body is constructed of stainless steel or material which is non-corrosive or has a corrosion resistant coating.

Designed for frequent and easy dismount and disinfection with hospital-grade products.

The design of the device should be such that it is easy to clean/disinfect by avoiding difficult-to-reach areas and having smooth surfaces.

---

The device has a sturdy design and shock resistant casing.



**Laboratory Instrument India**