



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Face mask, venturi, adult, single use

**Product Code :**

LBNY-0017-2200019



www.laboratoryinstrumentindia.com

**Description :**

Allows for precise oxygen delivery by utilizing different sized ports to change FiO<sub>2</sub> (24 – 60%). Face mask to provide oxygen at flow rates up to 12 LPM. Also known as air-entrainment masks able to provide total inspiratory flow at a specified FiO<sub>2</sub>.

**Technical Specification :**

Allows for precise oxygen delivery, with specific concentrations ranging from 24-60%.

The kit of the mask includes the tubing, humidity cup and multiple jets, which are color-coded and indicate the percentage of oxygen.

Mask is soft, transparent, with well-fitting mould.

The nose clip is adjustable.

Material: Mask and tubing PVC.

The tubing is terminated with standard or universal hose end connector.

Individually packed.

Tubing compatible with standard oxygen connection tubing.

Tubing length: 2 meters

Single use.

Non-sterile.

Adult size.

**ENVIRONMENTAL CONDITIONS**

Storage conditions: 0 to 60 deg. C, 10 to 90% relative humidity

Operating conditions: 5 to 40 deg. C, 10 to 85% relative humidity

---

**WEIGHT/VOLUME**

Volume: 0.001008m<sup>3</sup> (each, packaged)

Weight: 0.12kg (each, packaged).



**Laboratory Instrument India**