



Email : sales@laboratoryinstrumentindia.com

Product Name :

Face mask, venturi, child, single use

Product Code :

LBNY-0017-2200020



www.laboratoryinstrumentindia.com

Description :

Face mask to provide oxygen at flow rates up to 12 LPM. Allows for precise oxygen delivery by utilizing different sized ports to change FiO₂ (24 – 60%). Also known as air-entrainment masks able to provide total inspiratory flow at a specified FiO₂.

Allows for precise oxygen delivery, with specific concentrations ranging from 24-60%.

The kit of the mask includes the tubing, humidity cup and multiple jets, which are color-coded and indicate the percentage of oxygen.

Mask is soft, transparent, with well-fitting mould.

The nose clip is adjustable.

The tubing is terminated with standard or universal hose end connector.

Individually packed.

Tubing compatible with standard oxygen connection tubing.

Tubing length: 2 meters

Non-sterile.

Child size.

Single use.

Technical Specification :

ENVIRONMENTAL CONDITIONS

Storage conditions: 0 to 60 deg. C, 10 to 85% relative humidity

Operating conditions: 5 to 40 deg.C, 10 to 85% relative humidity

WEIGHT/VOLUME

Volume: 1.008dm³ (each, packaged)

Weight: 0.1kg (each, packaged).



Laboratory Instrument India