



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
High Speed Cutting M/C 180

**Product Code :**  
LBNY-0004-820022



**Description :**

High Speed Cutting M/C 180

**Technical Specification :**

High Speed Cutting M/C 180  
Saw blade size: 14 inch 66T  
No-load rotation speed: 1,300 RPM  
Motor: 15A 4HP  
Capacity (round): 5-3/16 inches  
Capacity (rectangular) : 4 - 1 /2 inches\* 6-1/2 inches  
Accessory  
- 10, 14 inch saw blades  
- 10 diamond blades  
- 1 hex key1 vertical clamp



**Laboratory Instrument India**

