



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Humidifier bottle, nonheated, reusable

**Product Code :**

LBNY-0017-2200024



www.laboratoryinstrumentindia.com

**Description :**

A sealed bottle of water inserted into a breathing circuit to add moisture to the breathing gases for administration to a patient. This type of humidifier does not heat the gas.  
Graduated, transparent humidification bottle, shatter resistant.  
Graduation shall show minimum and maximum water level.  
Detachable metal or rigid durable polymer cap with gas connectors.  
Non-heated, reusable humidifier for oxygen therapy and ventilation/anaesthesia inspiratory lines.  
Bubble-through humidification system.  
DISS connectors for gas inlet and outlet  
All available adaptors for gas connection to be described in bid response  
Flow rate capacity up to 15 L/min.  
Safety pop-off valve incorporated.  
Autoclavable and/or disinfectable with hospital grade detergents.  
Humidification chamber working volume at least 150 ml, not greater than 500 ml.

**Technical Specification :**

**ENVIRONMENTAL CONDITIONS**

Storage conditions temperature: -20 °C to 60°C / RH: 15% - 90%

Operating conditions temperature: 5 to 45°C / RH: 15 to 90%

**WEIGHT AND VOLUME:**

Volume: 0.627 dm<sup>3</sup> (per piece, packaged)

Weight: 0.146 kg (per piece, packaged).



**Laboratory Instrument India**