



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Sawing M/C 300

**Product Code :**  
LBNY-0004-820023



**Description :**

Sawing M/C 300

**Technical Specification :**

Sawing M/C 300  
Cutting depth: 41.5 mm  
Cable length: 3m  
Rotation speed (n): 3500 min<sup>-1</sup>  
Rated power (P<sub>n</sub>): 600W  
Cutting disc diameter: 120 mm  
Cutting disc dimensions: 120\*10\*2 mm -24T  
L<sub>p</sub>A sound pressure level: 94 dB(A)  
L<sub>w</sub>A sound power level: 105 dB(A)  
Rated AC voltage (One): 220-240 V-  
Rated frequency (ft): 50Hz  
Saw blade 10Ea



**Laboratory Instrument India**

