



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sawing M/C 300

Product Code :
LBNY-0004-820023



Description :

Sawing M/C 300

Technical Specification :

Sawing M/C 300
Cutting depth: 41.5 mm
Cable length: 3m
Rotation speed (n): 3500 min⁻¹
Rated power (P_n): 600W
Cutting disc diameter: 120 mm
Cutting disc dimensions: 120*10*2 mm -24T
L_pA sound pressure level: 94 dB(A)
L_wA sound power level: 105 dB(A)
Rated AC voltage (One): 220-240 V-
Rated frequency (ft): 50Hz
Saw blade 10Ea



Laboratory Instrument India

