



Email : sales@laboratoryinstrumentindia.com

Product Name :
Foam wraps, for YI sensor

Product Code :
LBNY-0017-2200037



www.laboratoryinstrumentindia.com

Description :

STORAGE INSTRUCTIONS

Storage Humidity: 10% to 95% relative humidity (non-condensing)

Storage Temperature: -40°C to +70°C

Technical Specification :

WEIGHT AND VOLUME

Volume: 0.49 dm³

Weight: 0.07 kg



Laboratory Instrument India