



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Foam wraps, for YI sensor

**Product Code :**  
LBNY-0017-2200037



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

**STORAGE INSTRUCTIONS**

Storage Humidity: 10% to 95% relative humidity (non-condensing)

Storage Temperature: -40°C to +70°C

**Technical Specification :**

**WEIGHT AND VOLUME**

Volume: 0.49 dm<sup>3</sup>

Weight: 0.07 kg



**Laboratory Instrument India**