



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Sensor, reusable, Rad-G, child/adult

**Product Code :**

LBNY-0017-2200038



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The Rad-G sensor is only compatible with the Rad-G device.

Prior to using this sensor the user should read and understand the Operator's Manual for the Rad-G device and the Directions for Use included with the reusable sensor.

The reusable sensors are indicated for noninvasive spot-checking of functional oxygen saturation of arterial hemoglobin (SpO<sub>2</sub>) and pulse rate (measured by an SpO<sub>2</sub> sensor) for use with adult, pediatric, and infant patients during both no motion and motion conditions, and for patients who are well or poorly perfused in hospitals, hospital-type facilities, mobile, and home environments.

**Technical Specification :**

**STORAGE INSTRUCTIONS**

Storage Humidity: 10% to 95% relative humidity (non-condensing)

Storage Temperature: -40°C to +70°C.

**WEIGHT AND VOLUME**

Volume: 0.49 dm<sup>3</sup>

Weight: 0.09 kg.

---

# Laboratory Instrument India

