



Email : sales@laboratoryinstrumentindia.com

Product Name :

Sensor, reusable, Rad-G, child/adult

Product Code :

LBNY-0017-2200038



www.laboratoryinstrumentindia.com

Description :

The Rad-G sensor is only compatible with the Rad-G device.

Prior to using this sensor the user should read and understand the Operator's Manual for the Rad-G device and the Directions for Use included with the reusable sensor.

The reusable sensors are indicated for noninvasive spot-checking of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (measured by an SpO₂ sensor) for use with adult, pediatric, and infant patients during both no motion and motion conditions, and for patients who are well or poorly perfused in hospitals, hospital-type facilities, mobile, and home environments.

Technical Specification :

STORAGE INSTRUCTIONS

Storage Humidity: 10% to 95% relative humidity (non-condensing)

Storage Temperature: -40°C to +70°C.

WEIGHT AND VOLUME

Volume: 0.49 dm³

Weight: 0.09 kg.

Laboratory Instrument India

