



Email : sales@laboratoryinstrumentindia.com

Product Name :

Sensor, reusable, Rad-G, neonate, YI set

Product Code :

LBNY-0017-2200039



www.laboratoryinstrumentindia.com

Description :

The Reusable Sensor is indicated for the continuous noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (measured by an SpO₂ sensor) for use with adult, pediatric, infant, and neonatal patients during both no motion and motion conditions and for patients who are well or poorly perfused in hospitals, hospital-type facilities, mobile, and home environments.

Technical Specification :

STORAGE INSTRUCTIONS

Storage Humidity: 10% to 95% relative humidity (non-condensing)

Storage Temperature: -40°C to +70°C.

WEIGHT AND VOLUME

Volume: 0.86 dm³

Weight: 0.12 kg.



Laboratory Instrument India

