



Email : sales@laboratoryinstrumentindia.com

Product Name :
Double End Grinder 450

Product Code :
LBNY-0004-820024



Description :

Double End Grinder 450

Technical Specification :

Double end grinder 450
8" (205mm) bench grinder
Amps: 4.2 amps
Peak horsepower: 3/4 hp
Spindle diameter: 5/8 in
Wheel diameter: 8"
Motor runs at 3,600 rpm for high-speed material removal
Bench grinder, crimp wire, 3 pcs for 8"
Wire wheel, 8", crimped, 5/8" arbor, wide face: 3 pcs.



Laboratory Instrument India