



Email : sales@laboratoryinstrumentindia.com

Product Name :
Resuscitator,hand-oper.,adult,set

Product Code :
LBNY-0017-2200045



www.laboratoryinstrumentindia.com

Description :

Hand-operated in a horizontal position.
Portable, self-inflating, reusable.
Ventilation can be done with ambient air or with oxygen.
The resuscitator has an inlet valve with a nipple for oxygen tubing.
The resuscitator is equipped with a non-rebreather valve.
An oxygen reservoir bag should be included with a capacity between 2,500-2,600ml.
Resuscitator can be fully disassembled, is easy to clean, disinfect, and reassemble.
All parts manufactured from high-strength, long-life materials that require no special maintenance or storage conditions.
The resuscitator, all its components and accessories are free of natural rubber latex.
The resuscitator is autoclavable between 121°C and 134°C.
The capacity of the self-inflating ventilation bag is 1,600ml.

Technical Specification :

ENVIRONMENTAL CONDITIONS

Storage conditions: -30°C - 50°C / 40% – 90% RH
Operating conditions: -18°C - 45°C / 15% – 90% RH

WEIGHT AND VOLUME

Volume: 6.86 dm³
Weight: 1.00 kg.



Laboratory Instrument India