



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Resuscitator,hand-oper.,adult,set

**Product Code :**  
LBNY-0017-2200045



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

#### **Description :**

Hand-operated in a horizontal position.  
Portable, self-inflating, reusable.  
Ventilation can be done with ambient air or with oxygen.  
The resuscitator has an inlet valve with a nipple for oxygen tubing.  
The resuscitator is equipped with a non-rebreather valve.  
An oxygen reservoir bag should be included with a capacity between 2,500-2,600ml.  
Resuscitator can be fully disassembled, is easy to clean, disinfect, and reassemble.  
All parts manufactured from high-strength, long-life materials that require no special maintenance or storage conditions.  
The resuscitator, all its components and accessories are free of natural rubber latex.  
The resuscitator is autoclavable between 121°C and 134°C.  
The capacity of the self-inflating ventilation bag is 1,600ml.

#### **Technical Specification :**

##### **ENVIRONMENTAL CONDITIONS**

Storage conditions: -30°C - 50°C / 40% – 90% RH  
Operating conditions: -18°C - 45°C / 15% – 90% RH

##### **WEIGHT AND VOLUME**

Volume: 6.86 dm<sup>3</sup>  
Weight: 1.00 kg.



**Laboratory Instrument India**