



Email : sales@laboratoryinstrumentindia.com

Product Name :
Nasal Catheter, 8 Fr

Product Code :
LBNY-0017-2200046



www.laboratoryinstrumentindia.com

Description :

Flexible nasal catheter with hole at distal end.
Oxygen and air/oxygen mixture compatibility.
Transparent tubing, suitable for low pressure gas supply.
Consists of connector, tubing and sponge at distal end.
Anti-kink tubing, non-permanent deformation if kinked or bent too tight.
Rubber or soft plastic tubing and catheter, semi-rigid and allowing freedom of movement, PVC certified for medical use, hardness > 60 Shore A.
8 Fr tubing, catheter length 40 cm.

Technical Specification :

ENVIRONMENTAL CONDITIONS

Storage conditions (packaged): Store under cool & dry conditions. Do not expose to temperatures above 49°C.
Operating conditions: 5 to 45°C, 15 to 90% relative humidity, non-condensing.

WEIGHT/VOLUME

Volume: 0.132 dm³ (each, packaged)
Weight: 8.00 gr (each, packaged).

Laboratory Instrument India

