



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Nasal Catheter, 8 Fr

**Product Code :**  
LBNY-0017-2200046



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Flexible nasal catheter with hole at distal end.  
Oxygen and air/oxygen mixture compatibility.  
Transparent tubing, suitable for low pressure gas supply.  
Consists of connector, tubing and sponge at distal end.  
Anti-kink tubing, non-permanent deformation if kinked or bent too tight.  
Rubber or soft plastic tubing and catheter, semi-rigid and allowing freedom of movement, PVC certified for medical use, hardness > 60 Shore A.  
8 Fr tubing, catheter length 40 cm.

**Technical Specification :**

**ENVIRONMENTAL CONDITIONS**

Storage conditions (packaged): Store under cool & dry conditions. Do not expose to temperatures above 49°C.  
Operating conditions: 5 to 45°C, 15 to 90% relative humidity, non-condensing.

**WEIGHT/VOLUME**

Volume: 0.132 dm<sup>3</sup> (each, packaged)  
Weight: 8.00 gr (each, packaged).

---

# Laboratory Instrument India

