



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Tool Box 250

**Product Code :**  
LBNY-0004-820025



**Description :**

Tool Box 250

**Technical Specification :**

Tool Box 250  
3-tier Rolling Mobile Tool Storage Box  
1st stage (top)  
IP65 waterseal  
Metal wire front latches  
Toughsystem 1 1mouldlies  
2nd stage (middle)  
Ip65 water seal  
3rd tier (bottom)  
Detachable Handel  
8" pu durable Wheels  
Ip65 waterseal  
Name panel insert

---

# Laboratory Instrument India

