



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tool Box 250

Product Code :
LBNY-0004-820025



Description :

Tool Box 250

Technical Specification :

Tool Box 250
3-tier Rolling Mobile Tool Storage Box
1st stage (top)
IP65 waterseal
Metal wire front latches
Toughsystem 1 1moudlies
2nd stage (middle)
Ip65 water seal
3rd tier (bottom)
Detachable Handel
8" pu durable Wheels
Ip65 waterseal
Name panel insert

Laboratory Instrument India

