



Email : sales@laboratoryinstrumentindia.com

Product Name :
Drum,sterilizing,165mm diam

Product Code :
LBNY-0018-220003



www.laboratoryinstrumentindia.com

Description :

Cylindrical container used to sterilise dressing materials (gauze compress or cotton etc.) in a steam sterilizer (autoclave), and to keep them as "sterile" dressing materials for medical activities (i.e. dressing, injection etc.)
Drum, sterilizer.

Air vents system to allow steam to circulate freely during the sterilisation cycle.

Lateral air vents system is preferable to top and bottom air vents.

Closing lid with a clip lock.

Material: Austenitic stainless steel, smooth surface.

External diameter: approx. 150 to 165 mm.

Height: approx. 100 to 120 mm.

Thickness: approx. 0.6 to 0.7 mm.

Carrying handle.

Austenitic stainless steel composition: approx. 8-10 % nickel, 18-20 % chromium.

Reusable.

Technical Specification :

Environmental conditions:

Humidity: 40% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.

Laboratory Instrument India

