



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Drum,sterilizing,260mm diam

**Product Code :**  
LBNY-0018-220004



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Drum, sterilizer.

Air vents system to allow steam to circulate freely during the sterilisation cycle.

Lateral air vents system is preferable to top and bottom air vents.

Material: Austenitic stainless steel, smooth surface.

Closing lid with a clip lock.

Carrying handle.

Height: approx. 160 to 170 mm.

Thickness: approx. 0.6 to 0.7 mm.

External diameter: approx. 240 to 260 mm.

Austenitic stainless steel composition: approx. 8-10 % nickel, 18-20 % chromium.

Reusable.

**Technical Specification :**

Environmental conditions:

Humidity: 40% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.

---

# Laboratory Instrument India

