



Email : sales@laboratoryinstrumentindia.com

Product Name :
Drum,sterilizing,340mm diam

Product Code :
LBNY-0018-220006



www.laboratoryinstrumentindia.com

Description :

Drum, sterilizer.

Material: Austenitic stainless steel, smooth surface.

Austenitic stainless steel composition: approx. 8-10 % nickel, 18-20 % chromium.

Reusable.

Air vents system to allow steam to circulate freely during the sterilisation cycle.

Lateral air vents system is preferable to top and bottom air vents.

Closing lid with a clip lock.

Carrying handle.

External diameter: approx. 340 to 350 mm.

Height: approx. 180 to 200 mm.

Thickness: approx. 0.6 to 0.7 mm.

Technical Specification :

Environmental conditions:

Humidity: 40% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.

Laboratory Instrument India

