



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vertical Drill

Product Code :
LBNY-0004-820026



Description :

Vertical Drill

Technical Specification :

Vertical Drill
Head is honed and oiled parts are cylindrical and flat.
With advanced keyless chuck
Belt cover with a thickness of 1mm or more
Power:220V / singlephase
Drilling ability: 19mm
Chuck Specification: MT2
Swing; 180mm
Curling diameter: 80mm
Distance between main shaft end and base up to 650mm
Up and down movement distance: 85mm
Drill bit: 10Ea



Laboratory Instrument India

