



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Masking tape,for ster. pack

**Product Code :**  
LBNY-0018-220009



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Used to close paper crepe packs for steam sterilization. Easy released pressure sensitive adhesive, easy to tear paper, easy to remove.

Pressure sensitive adhesive

Paper-based adhesive tape

Leaves no residue

Does not damage the surface to which it is applied.

Plain (without sterilization indicators)

Maintains integrity and sticks to surface during steam sterilization cycle and drying temperatures.

Size: width 17-21mm x length 48-52m.

**Technical Specification :**

Weigh and volume:

Estimated volume: 0.355CDM

Estimated weight: 0.080Kg



**Laboratory Instrument India**

