



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Stove,kerosene,single-burner,pressure

**Product Code :**

LBNY-0018-2200010



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Stove, kerosene, single-burner, pressure

Kerosene stove designed to provide safe and sufficient heat for medical sterilization applications.

Hand-operated plunger pump allows pressurisation of fuel tank and control of flame intensity

Pressure, max: 1.3 bar

Comprises of a fuel tank, vapour burner and pot holder

Single burner with top plate

Fit with overpressure safety valve

Material: Heat and rustproof, steel or brass

Capacity kerosene fuel tank: 1.2 to 2.5 L

Diameter flame holder: 5-7 cm

Fuel consumption: 0.4-0.6 L/h

Height: 20-29 cm

**Technical Specification :**

Weigh and volume:

Estimated volume: 28.196CDM

Estimated weight: 2.855Kg.

---

# Laboratory Instrument India

